

ICX-ALL-IN-ONE[®] DRILLS

MORE THAN
400 IMPLANTS
WITH **1 DRILL!**

PATENTED
WORLD INNOVATION!



DISCOVER THE ICX-DIFFERENCE.


medentis
medical

ICX-ALL-IN-ONE[®] DRILLS

ICX-ALL-IN-ONE-Drills

The ICX-ALL-IN-ONE DRILLS gradually expand the diameter during the surgical procedure and are designed to be used with standard surgical equipment to precisely cut the bone with clockwise rotation (800 rpm/ pre-drill and 400 rpm/normal drill).

Clockwise,
cutting direction



ICX-SHOP!

Cutting mode



www.medentis.de

medentis
medical

ICX-ALL-IN-ONE[®] DRILLS



*The World
Innovation*



ICX-SHOP!

www.medentis.de

medentis
medical


Ø 2.9mm D2/D3/D4

Drill protocol for very soft/
soft/medium hard bone.

This procedure is identical
for all lengths.

D2/D3/D4 drill,
1 ring

Cortical drill



The drilling protocol suggested by medentis medical is only a recommendation.

In general, the rule for all drills, for all implant diameters and in all bone classes is, that the prescribed drill sequence can vary. If the bone is softer, fewer steps are recommended. The harder the bone, the more intermediate steps are recommended.

ICX-ALL-IN-ONE
DRILLS

Ø 2.9mm D1

Drill protocol for hard bone. Please use only drills with 3 black rings! This procedure is identical for all lengths.

Cortical drill



D1 drill
3 rings,



medentis
medical

www.medentis.de

The drilling protocol suggested by medentis medical is only a recommendation.

In general, the rule for all drills, for all implant diameters and in all bone classes is, that the prescribed drill sequence can vary. If the bone is softer, fewer steps are recommended. The harder the bone, the more intermediate steps are recommended.


Ø 3.3mm D2/D3/D4

Drill protocol for very soft/
soft/medium hard bone.

This procedure is identical
for all lengths.

D2/D3/D4 drill, 
1 ring

Cortical drill



The drilling protocol suggested by medentis medical is only a recommendation.

In general, the rule for all drills, for all implant diameters and in all bone classes is, that the prescribed drill sequence can vary. If the bone is softer, fewer steps are recommended. The harder the bone, the more intermediate steps are recommended.


ICX-ALL-IN-ONE[®]
DRILLS

Ø 3.3mm D1

Drill protocol for hard bone. Please only use drills with 3 pink rings here! This procedure is identical for all lengths.

Cortical drill

D1 drill
3 rings,



medentis
medical

www.medentis.de

The drilling protocol suggested by medentis medical is only a recommendation.

In general, the rule for all drills, for all implant diameters and in all bone classes is, that the prescribed drill sequence can vary. If the bone is softer, fewer steps are recommended. The harder the bone, the more intermediate steps are recommended.

Ø 3.45mm D2/D3/D4

Drill protocol for very soft/
soft/medium hard bone.

This procedure is identical
for all lengths.


The drilling protocol suggested by medentis medical is only a
recommendation.

In general, the rule for all drills, for all implant diameters and in
all bone classes is, that the prescribed drill sequence can vary. If
the bone is softer, fewer steps are recommended. The harder the
bone, the more intermediate steps are recommended.

ICX-ALL-IN-ONE[®]
DRILLS

D2/D3/D4 drill, ←

Cortical drill



Ø 3.45mm D1

Drill protocol for hard bone. Please only use drills with 3 yellow rings here! This procedure is identical for all lengths.

Cortical drill



D1 drill
3 rings,



medentis
medical

www.medentis.de

The drilling protocol suggested by medentis medical is only a recommendation.

In general, the rule for all drills, for all implant diameters and in all bone classes is, that the prescribed drill sequence can vary. If the bone is softer, fewer steps are recommended. The harder the bone, the more intermediate steps are recommended.


Ø 3.75mm D2/D3/D4

Drill protocol for very soft/
soft/medium hard bone.

This procedure is identical
for all lengths.

D2/D3/D4 drill, 1 ring ←

Cortical drill



The drilling protocol suggested by medentis medical is only a recommendation.

In general, the rule for all drills, for all implant diameters and in all bone classes is, that the prescribed drill sequence can vary. If the bone is softer, fewer steps are recommended. The harder the bone, the more intermediate steps are recommended.

ICX-ALL-IN-ONE[®]
BOHRER


Ø 3.75mm D1

Drill protocol for hard bone. Please only use drills with 3 red rings here! This procedure is identical for all lengths.

Cortical drill

D1 drill
3 rings, red

Optional
to use



medentis
medical

www.medentis.de

The drilling protocol suggested by medentis medical is only a recommendation.

In general, the rule for all drills, for all implant diameters and in all bone classes is, that the prescribed drill sequence can vary. If the bone is softer, fewer steps are recommended. The harder the bone, the more intermediate steps are recommended.


Ø 4.1mm D2/D3/D4

Drill protocol for very soft/
soft/medium hard bone.

This procedure is identical
for all lengths.

Cortical drill

D2/D3/D4 drill, 1 ring ←



The drilling protocol suggested by medentis medical is only a recommendation.

In general, the rule for all drills, for all implant diameters and in all bone classes is, that the prescribed drill sequence can vary. If the bone is softer, fewer steps are recommended. The harder the bone, the more intermediate steps are recommended.

ICX-ALL-IN-ONE
DRILLS

Ø 4.1mm D1

Drill protocol for hard bone. Please only use drills with 3 green rings here! This procedure is identical for all lengths.

Cortical drill

D1 drill
3 rings, green

Optional
to use

medentis
medical

www.medentis.de

The drilling protocol suggested by medentis medical is only a recommendation.

In general, the rule for all drills, for all implant diameters and in all bone classes is, that the prescribed drill sequence can vary. If the bone is softer, fewer steps are recommended. The harder the bone, the more intermediate steps are recommended.


Ø 4.8mm D2/D3/D4

Drill protocol for very soft/
soft/medium hard bone.

This procedure is identical
for all lengths.

D2/D3/D4 drill, ←
1 ring

Cortical drill



The drilling protocol suggested by medentis medical is only a recommendation.

In general, the rule for all drills, for all implant diameters and in all bone classes is, that the prescribed drill sequence can vary. If the bone is softer, fewer steps are recommended. The harder the bone, the more intermediate steps are recommended.

ICX-ALL-IN-ONE[®]
DRILLS

Ø 4.8mm D1

Drill protocol for hard bone. Please only use drills with 3 blue rings here! This procedure is identical for all lengths.

Cortical drill

D1 drill
3 rings, green

D1 drill
3 rings, blue

Optional
to use



medentis
medical

www.medentis.de

The drilling protocol suggested by medentis medical is only a recommendation.

In general, the rule for all drills, for all implant diameters and in all bone classes is, that the prescribed drill sequence can vary. If the bone is softer, fewer steps are recommended. The harder the bone, the more intermediate steps are recommended.

ICX-ALL-IN-ONE[®]


SURGERY BOX

*The World
Innovation*



VORSPRUNG DURCH INNOVATION.

medentis
medical



ICX-ALL-IN-ONE[®]

SURGERY BOX



ICX-SHOP!

ICX-ALL-IN-ONE SURGERY BOX | TABLE OF CONTENTS

Pos.-No.		Article-No.	Name	Discription
1	2	AIO-014-200000	All in one Pre-drill	Ø 2,0
2	1	AIO-014-480049	All in one drill	Ø 4,8 x 4,9
3	1	AIO-014-345065	All in one Drill	Ø 3,45 x 6,5
4	1	AIO-014-480065	All in one Drill	Ø 4,8 x 6,5
5	1	AIO-014-290080	All in one Drill	Ø 2,9 x 8
6	1	AIO-014-330080	All in one Drill	Ø 3,3 x 8
7	1	AIO-014-345080	All in one Drill	Ø 3,45 x 8
8	1	AIO-014-375080	All in one Drill	Ø 3,75 x 8
9	1	AIO-014-410080	All in one Drill	Ø 4,1 x 8
10	1	AIO-014-480080	All in one Drill	Ø 4,8 x 8
11	1	AIO-014-290100	All in one Drill	Ø 2,9 x 10
12	1	AIO-014-330100	All in one Drill	Ø 3,3 x 10
13	1	AIO-014-345100	All in one Drill	Ø 3,45 x 10
14	1	AIO-014-375100	All in one Drill	Ø 3,75 x 10
15	1	AIO-014-410100	All in one Drill	Ø 4,1 x 10
16	1	AIO-014-480100	All in one Drill	Ø 4,8 x 10
17	1	AIO-014-290125	All in one Drill	Ø 2,9 x 12,5
18	1	AIO-014-330125	All in one Drill	Ø 3,3 x 12,5
19	1	AIO-014-345125	All in one Drill	Ø 3,45 x 12,5
20	1	AIO-014-375125	All in one Drill	Ø 3,75 x 12,5
21	1	AIO-014-410125	All in one Drill	Ø 4,1 x 12,5
22	1	AIO-014-480125	All in one Drill	Ø 4,8 x 12,5
23	1	AIO-014-375150	All in one Drill	Ø 3,75 x 15
24	1	AIO-014-410150	All in one Drill	Ø 4,1 x 15
25	1	AIO-014-480150	All in one Drill	Ø 4,8 x 15
26	1	AIO-014-007290	All in one Drill	Ø 2,9 D1
27	1	AIO-014-007330	All in one Drill	Ø 3,3 D1
28	1	AIO-014-007345	All in one Drill	Ø 3,45, D1
29	1	AIO-014-007375	All in one Drill	Ø 3,75, D1
30	1	AIO-014-007410	All in one Drill	Ø 4,1, D1
31	1	AIO-014-007480	All in one Drill	Ø 4,8, D1

32	1	AIO-014-004330	All in one Thread cutter	Ø 3,3
33	1	AIO-014-004375	All in one Thread cutter	Ø 3,75
34	1	AIO-014-004410	All in one Thread cutter	Ø 4,1
35	1	AIO-014-004480	All in one Thread cutter	Ø 4,8
36	4	C-014-350410	Parallel pin	D=3.5/4.1mm
37	4	C-014-230315	Parallel pin	D=2.3/3.15mm
38	1	960004	Drill extension	
39	1	C-014-006002	Bone cutter	
40	1	C-014-006003	ICX 3.3 Bone cutter	rosa / rose
41	1	C-014-000003	Crestotom Drill	
42	1	C-015-100007	ICX t-bona insertion tool ISO shaft	ICX T-bona screw
43	1	C-020-000001	Multi insertion tool ISO shaft	L=24mm
44	1	960002	Screwdriver ISO-shaft	hex 0.9mm L=20mm
45	1	C-015-100006	ICX-Mini insertion tool ISO-shaft	
46	1	C-015-100033	Screwdriver ISO-shaft	hex 1.2mm L=25mm
47	1	N-015-100002	ICX-3.3 insertion tool ISO-shaft	L=29mm
48	1	950099	Screwdriver ISO-shaft	hex 1.4mm short
49	1	C-015-100002	ICX insertion tool ISO-shaft	L=29mm
50	2	960007	Universal ratchet adapter	
51	1	C-015-100029	Screwdriver	hex 0.9mm
52	1	C-015-100030	Screwdriver	hex 1.2mm
53	1	C-015-100025	Screwdriver	hex 1.4mm mittel/middle
54	1	960001	Torque ratchet	
55	1	C-015-100009	ICX-Mini Counter tool	
56	1	C-015-100003	Depth measuring probe	
57	1	C-018-100000	ICX All in one Box	not filled

ICX-ALL-IN-ONE[®] DRILLS

MORE THAN
400 IMPLANTS
WITH **1 DRILL!**



ICX-SHOP!

PATENTED
WORLD INNOVATION!

DISCOVER THE ICX-DIFFERENCE.

medentis
medical

www.medentis.de